



# Where is your AED...?

In the event of an emergency...knowing where, and having access to your Automated External Defibrillator can be a matter of Life or Death.

*"Each year, more than 250,000 Americans die of sudden out-of-hospital cardiac arrest. Survival from sudden cardiac arrest in adults depends directly on the speed of defibrillation; every minute of delay in defibrillation reducing the chances of survival 7 to 10%".*

Virginia Department of Health—Office of Emergency Medical Services





The Series 1750 Wall Mount Enclosure was not simply designed to conveniently store your AED; it was designed as an integral part of your lifesaving Public Access Defibrillation Program.

[Click to order](#)

When sudden cardiac arrest occurs, time is of the essence. Trained emergency support personnel must be alerted to the situation. By simply opening the access door to the Series 1750 Wall Mount Enclosure, a 110 db., audible alarm will begin to sound, a frequency loud enough to alert



everyone in the area to the emergency situation. A strobe option is also available to visually pinpoint the exact location of the AED needed for use.



If a central security office is located on the premises, a hard wire hook-up to the Series 1750 Wall Mount Enclosure is offered as an option and can be used to notify others outside the immediate area of emergency situation.

Since restricting access can prove deadly in case of an emergency, the Series 1750 Wall Mount Enclosure is equipped with a simple non-locking push button door latch. However, the same lifesaving warning devices can also serve to protect your AED from tampering and theft.

**AED Location Placard** - Heavy gauge ridged plastic placard - clear signage protruding about five inches from the wall indicating your AED location.



[Click to order](#)